

NARISHIGE

GDC-1/1.2/1.5

Glass Capillary with Filament

Box of 200 (5 each to a vacuum packed bag).

These are glass capillaries with a filament inside. The filament leads the liquid up to the capillary tip. They are pre-washed in an ultrasonic washing machine, sterilized with ultraviolet ray, and vacuum-packed.



NARISHIGE

GDC-1/1.2/1.5

Glass Capillary with Filament

Specification

Type	<u>G-1</u>	<u>G-1.2</u>	<u>G-1.5</u>	<u>G-2</u>	<u>G-3</u>	<u>GD-1</u>	<u>GD-1.2</u>	<u>GD-1.5</u>	<u>GD-2</u>
Outer diameter (mm)	1	1.2	1.5	2	3	1	1.2	1.5	2
Inner diameter (approx. mm)	0.6	0.75	0.9	0.8	1.8	0.6	0.75	0.9	1.2
Length (mm)	90*	90*	90*	90*	90	90*	90*	90*	90*
Contents (pcs)	500	500	500	200	200	500	500	500	200
Material	Borosilicate								

Type	<u>GC-1</u>	<u>GC-1.2</u>	<u>GC-1.5</u>	<u>GDC-1</u>	<u>GDC-1.2</u>	<u>GDC-1.5</u>	<u>G-100</u>	<u>G-100L</u>	<u>G-1000</u>
Outer diameter (mm)	1	1.2	1.5	1	1.2	1.5	1	1	1
Inner diameter (approx. mm)	0.6	0.75	0.9	0.6	0.75	0.9	0.75	0.75	-
Length (mm)	90	90	90	90	90	90	90*	100*	90*
Contents (pcs)	200	200	200	200	200	200	500	500	200
Material	Borosilicate								

* A custom length between 70mm and 180mm is also available by special order.

Specifications are subject to change without notice.

Copyright © NARISHIGE Group. All Rights Reserved. 1999-2020