



# MAG-1

## Simple Head Holder Plate (for Mice)

### For observations of a mouse under a microscope.

After bonding a chamber onto the head of a mouse on the SR-9M-HT Stereotaxic Instruments, the MAG-1 can be transferred and placed under a microscope or other experimental instrument for performing observations/experiments. Designed to be compact in size, the MAG-1 can be placed under a two-photon microscope or stereomicroscope.



\* Attachable chambers: CF-10, CP, CP-1, CP-2


## Specification

<b>Accessory included</b>	<u>CF-10</u> Chamber Frame (1pcs)
<b>Dimensions/Weight</b>	W96 × D84 × H42mm, Approx. 102g

\* The chamber frame is customizable to your own application. Please contact us for details.

\* The CFR-1 / CFR-2 Chamber Frame (for rats) can NOT be attached.

## System Diagram



● refers to the model itself with standard attachments

○ refers to attachable main accessories

○ refers to types of mounting adaptors

Specifications are subject to change without notice.

Copyright © NARISHIGE Group. All Rights Reserved. 1999-2019